

Comparison of RNP-mediated editing by Type V Cpf1 variants across multiple cell types

John Zuris, Ramya Viswanathan, Katherine Loveluck, Elise Keston-Smith, Eric Tillotson, Kiran Gogi, Jack Heath, Aditi Chalise, Vidya Dhanapal, Tongyao Wang, Melissa Chin, Justin Fang, Chris Borges, Luis Barrera, Frank Buquicchio, Carrie Margulies, Anne Bothmer, Cecilia Cotta, Georgia Giannoukos, Grant Welstead, Fred Harbinski, Jennifer Gori, Chris Wilson, Vic Myer, Hari Jayaram

Abstract

The CRISPR-Cpf1 system offers several potential advantages over other nucleases for *ex vivo* genome editing therapies, including a smaller single crRNA that can be readily synthesized, the ability to target T- and C-rich PAMs with the WT and RR variants, respectively, and lastly a 5'-staggered cut which may lead to different repair outcomes (1).


We have optimized several Type V Cpf1 variant ribonucleoproteins which are the preferred delivery mode for *ex vivo* gene-editing therapeutics. Comparing their cellular potency with SpCas9 we show that multiple Cpf1 variants show robust editing activity at multiple sites in cells.

In addition, we show that Cpf1 orthologs such as Lb2Cpf1 and FnCpf1 can lead to robust editing in cells when delivered as RNPs.


Lastly, we demonstrate efficient editing with several Cpf1 variants delivered as RNPs in T cells and AsCpf1 delivered as an RNP in adult hematopoietic stem cells (HSCs).

These findings underscore the promise of RNP delivery for Cpf1 nucleases which have desirable properties for *ex vivo* genome editing therapeutics.


Robust pipeline for production and evaluation of our RNPs




Efficient editing with Cpf1 orthologs delivered as RNPs




Engineered Cpf1 variants expand PAM targeting space



Efficient editing in T cells at multiple loci with Cpf1 RNPs



Efficient editing in HSCs at multiple loci with Cpf1 RNPs



Cpf1 vs. Cas9

Cpf1	Cas9
<ul style="list-style-type: none"> WT TTTV, TTN PAMs plus TYCV/CCCC, TATV engineered PAM variants Single guide with 20-24 nt protospacer and 19-20 nt direct repeat (~40 nt) 5' staggered DNA cut with 4-5 nt overhangs 	<ul style="list-style-type: none"> WT NGG PAM plus NGAN, NGCG engineered PAM variants Separate crRNA and tracrRNA that can be linked into one molecule (~100 nt) Blunt DNA cut

References:

1. Zetsche et al. *Cell* 2015
2. Kleinstiver et al. *Nat Biotech* 2016
3. Gao et al. *Biorxiv* 2016