

# An in vitro method for benchmarking of CRISPR-associated endonucleases



Eric Tillotson<sup>1</sup>, Gregory Gotta<sup>1</sup>, Katherine Loveluck<sup>1</sup>, John Zuris<sup>1</sup>, Barrett Steinberg<sup>1</sup>, Tongyao Wang<sup>1</sup>, Ramya Viswanathan<sup>1</sup>, Elise Keston-Smith<sup>1</sup>, Joost Mostert<sup>1</sup>, Luis Barrera<sup>1</sup>, Christopher Wilson<sup>1</sup>, Vic Myer<sup>1</sup>, Hari Jayaram<sup>1</sup>

<sup>1</sup>Editas Medicine, 11 Hurley Street, Cambridge, MA. 02141


## Introduction

Direct delivery of Cas9 or Cpf1 ribonucleoprotein (RNP) complexes to cells by electroporation is emerging as an important delivery mode for ex vivo gene editing therapies. Understanding the dose/efficacy relationship of the delivered RNP is a key step in optimizing a therapeutic. An essential requirement in accurately establishing this relationship and determining RNP quality is measuring the *in vitro* activity of a given RNP complex in a standardized manner. Here we describe an *in vitro* assay for measuring the specific activity of CRISPR-associated nucleases on synthetic DNA substrates. We then apply this assay to benchmark the activity of several naturally occurring and engineered CRISPR nucleases. This standardized *in vitro* benchmarking provides critical input to help create the best RNP based therapy.

## Specific Activity Method


The specific activity assay substrate (a) is a 200 bp synthetic DNA amplicon comprised of genomic DNA target sequences and constant 5' and 3' primer handles. Selected target sites have a 3' NGGRR PAM with one of the following 5' PAMs 20 bp upstream: TTTV, TTN, NYC, or TATV, allowing us to compare activities of multiple nucleases using the same synthetic substrate. EC50 values are generated by a dose-response of RNP over 15 nM substrate for 1 hour at 37° C in 30 mM Hepes pH 7.5, 150 mM NaCl, and 2 mM MgCl<sub>2</sub>. RNP is added to DNA substrate using a Biomek FX<sup>9</sup> liquid handler. Nuclease activity is quenched by addition of Proteinase K. Fraction of cleaved substrate is measured using the Fragment Analyzer automated CE system (AATI). Analysis is performed with GraphPad Prism using four parameter dose-response.

## In vitro Benchmarking of CRISPR-associated Nucleases




a) Specific activity benchmarking was performed on the following CRISPR nucleases: SpCas9, SauCas9, High-fidelity (SpCas9-HF1), enhanced specificity (eSpCas9), AsCpf1, LbCpf1, FnCpf1, AsCpf1-RVR, and AsCpf1 RR. N = 2 b) 95% CI EC50 values indicate that the variance in activity between enzymes is similar and largely target specific. NGG PAM was used for SpCas9, SpCas9-HF1, and eCas9. NGGRRRT for SauCas9. TTTV for AsCpf1, LbCpf1, and FnCpf1. TATV for AsCpf1 RVR. NYC for AsCpf1 RR.


## Assay Characterization and Sensitivity



a) Assay sensitivity was determined by doping in of i. dead (orthogonal) gRNA and ii. dSpCas9 b) 95% confidence interval values indicate EC50 values are significantly different at 50% i. dead gRNA and ii. dSpCas9 c) Decrease in Amax suggests dSpCas9 binds to substrate DNA inhibiting SpCas9 activity.



## Off-target Analysis of Engineered Variants



a) SpCas9, SpCas9-HF1, and eSpCas9 show comparable on-target activity. Rationally engineered SpCas9-HF1 and eSpCas9 have almost no activity against off-targets *in vitro*.

## Conclusions

*In vitro* matched site benchmarking offers activity comparisons between CRISPR-associated nucleases. Potency measurements on matched substrates indicate no significant differences in activity across nucleases tested. Finally, engineered variants High-fidelity Cas9 and enhanced specificity Cas9 display target discrimination *in vitro*. Benchmarking these nucleases enables understanding the cellular dose-efficacy relationship and helps optimize an ex-vivo RNP therapeutic.

Disclaimer: This presentation contains forward-looking statements within the meaning of the "safe harbor" provisions of the Private Securities Litigation Reform Act of 1995. All statements, other than statements of historical facts, contained in this presentation, including statements regarding the Company's strategy, future operations, future financial position, future revenue, projected costs, prospects, plans, and objectives of management, are forward-looking statements. The words "anticipate," "believe," "continue," "could," "estimate," "expect," "intend," "may," "plan," "potential," "predict," "project," "should," "would," and similar expressions are intended to identify forward-looking statements, although not all forward-looking statements contain these identifying words. Forward-looking statements in this presentation include statements regarding the Company's goals of submitting an IND for the LCA10 program in 2017, entering an LCA10 natural history study in mid-2017, achieving preclinical proof-of-concept for additional programs and establishing a database, the Company may not actually achieve the plans, intentions, or expectations disclosed in these forward-looking statements, and you should not place undue reliance on these forward-looking statements. Actual results or events could differ materially from the plans, intentions and expectations disclosed in these forward-looking statements as a result of various factors, including: uncertainties inherent in the initiation and completion of preclinical studies and clinical trials and clinical development of the Company's product candidates; clinical trials whether interim results from a clinical trial will be predictive of the final results of the trial or the results of future trials; expectations for regulatory approvals to conduct trials or to market products; availability of funding sufficient for the Company's foreseeable and unforeseeable operating expenses and capital expenditure requirements; and other factors discussed in the "Risk Factors" section of the Company's most recent Quarterly Report on Form 10-Q, which is on file with the Securities and Exchange Commission, and in other filings that the Company may make with the Securities and Exchange Commission in the future. In addition, the forward-looking statements included in this presentation represent the Company's views as of the date of this presentation. The Company anticipates that subsequent events and developments will cause its views to change. However, while the Company may elect to update these forward-looking statements at some point in the future, it specifically disclaims any obligation to do so. These forward-looking statements should not be relied upon as representing the Company's views as of any date subsequent to the date of this presentation.